



Anti-Myc Tag Mouse Monoclonal Antibody (2D5), Cy3 Conjugated

Cat #: A02060CY3

Size: 100µl

Product Information

	Product Name: Anti-Myc Tag Mouse Monoclonal Antibody (2D5), Cy3 Conjugated		
	Applications: IF		Isotype: Mouse IgG
	Reactivity: All Species Expected		
REF	Catalog Number: A02060CY3	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C protected from light. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: Myc tag is a polypeptide protein tag derived from the c-myc gene product that can be added to a protein using recombinant DNA technology. It can be used for affinity chromatography, and then used to separate recombinant, overexpressed protein from wild type protein expressed by the host organism. It can also be used in the isolation of protein complexes with multiple subunits.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IF (1:1000).

Storage Buffer: Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in the dark.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.