


Human GDNF protein

Cat #: PRP1067

Size: 5µg/20µg/100µg/1mg

Product Information

	Product Name: Human GDNF protein		
REF	Catalog Number: PRP1067	LOT	Lot Number: Refer to product label
	Purity: > 95% as determined by SDS-PAGE		
	Storage: Store at -20°C		Preparation method: E.coli
	Shipping: The product is shipped at ambient temperature.		

Background: Glial cell line-derived neurotrophic factor (GDNF) is a glycosylated, disulfide-bonded homodimer that is a distantly related member of the transforming growth factor-β superfamily. It is a founding member of the GDNF family of ligands (GFL) and has been shown to interact with GFRA2 and GDNF family receptor alpha. Mature Rat and Human GDNF exhibit approximately 93% amino acid sequence identity and show considerable species cross-reactivity. Cells known to express GDNF include Sertoli cells, type 1 astrocytes, Schwann cells, neurons, pinealocytes and skeletal muscle cells. In embryonic midbrain cultures, recombinant Human GDNF promoted the survival and morphological differentiation of dopaminergic neurons and increased their high-affinity dopamine uptake.

Sequence: Amino acid sequence derived from Human GDNF protein (Ser78-Ile211) with 6×His tag at the C-terminus was expressed.

Protein length: The recombinant Human GDNF consists of 134 amino acids and predicts a molecular mass of 16 kDa.

Biological Activity: Measure by its ability to induce proliferation in SH-SY5Y cells. The ED₅₀ for this effect is <10 ng/mL.


Endotoxin: < 0.1 EU per µg of the protein as determined by the LAL method.

Formulation: The protein was lyophilized from 20 mM sodium citrate, 0.2 M NaCl, pH 3.5.

Storage Instructions: Lyophilized protein product should be stored desiccated below -18°C. Upon reconstitution, the protein should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (5% HSA, 10%FBS or 0.1%BSA). Please prevent freeze-thaw cycles.

Usage notes: Always centrifuge tubes before opening. It is recommended to reconstitute the lyophilized recombinant Human

GDNF protein to a concentration of 0.1-1 mg/mL in sterile ddH₂O, and keep at room temperature for at least 20 min to fully dissolve. Please avoid vortex vigorously.



Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.